

STAY CONNECTED WITH THE

REACT DATA LINK

WIRELESS DATA SHARING ACCESSORY



- Seamlessly share data with patient care team
- Easy remote access to:
 - View trend reports
 - Track therapy compliance
 - Facilitate informed treatment decisions
- Simple, reliable, and HIPAA Compliant
- GPS allows inventory traceability



Wireless Data
Sharing



HIPAA Compliant



GPS Location
Support



Easy USB Plug &
Play Functionality



COORDINATE CARE



WITH CONVENIENT COMMUNICATION

The React DataLink provides wireless support for patient data sharing, allowing clinicians to view trends, track patient compliance with prescribed therapy, and facilitate informed treatment decisions.



Simple, reliable, and HIPAA compliant, the React DataLink encrypts data to communicate over a cellular network. For customers with limited access to cellular networks, the React DataLink can store up to 30 days of data to share when it next connects to an appropriate network. In addition to data sharing, the React DataLink also provides GPS location support for connected devices.

The React DataLink is fully compatible with the hardware on any existing VOCSN product and can easily be attached in the field without tools. A simple software update will be required for any VOCSN operating with software older than version 5.00.09.

MULTI-VIEW CONNECT

Access to connected device data is managed on a dedicated, password protected web portal called Multi-View Connect. With unique accounts, qualified users can access appropriate device reports that can be customized for specific caregiver roles such as physician, respiratory therapist, or DME. A cellular network subscription is required to support access to Multi-View Connect.

To learn more, visit: venteclife.com.